

---

## limited 1 year warranty

---

Your Inspired by Drive branded product is warranted to be free of defects in materials and workmanship for 1 year from the date of purchase for the original consumer purchaser.

This device was built to exacting standards and carefully inspected prior to shipment. This Limited 1 Year Warranty is an expression of our confidence in the materials and workmanship of our products and our assurance to the consumer of years of dependable service.

This warranty does not cover device failure due to owner misuse or negligence, or normal wear and tear. The warranty does not extend to non-durable components, such as rubber accessories, casters, and grips, which are subject to normal wear and need periodic replacement.

If you have a question about your Inspired by Drive device or this warranty, please contact an authorized Inspired by Drive dealer.



2017 Medical Depot, Inc. All rights reserved.

Drive DeVilbiss Healthcare

99 Seaview Boulevard

Port Washington N.Y. 11050 USA

Made in China



## Kanga solid seat pan

Item #


KGA-1635-2G

KGA-1835-2G

KGA-2035-2G

[www.inspiredbydrive.com](http://www.inspiredbydrive.com)

## assembly



1. Remove knobs (A) from seat pan attaching brackets
2. Mount seat pan on seat rails with the cutout facing the rear of the chair.
3. Push seat pan onto seat rails making sure the mounting brackets (B) click in place.
4. Install connecting knobs (A)

